


*Cryo-Line*



**OPTIKON™**  
Man and Technology

# REF 131002

- Simple and automatic system for cryosurgery.
- Wide graphic display.
- User friendly control software.
- Probe temperature control.
- Count down function.
- High precision temperature reading.
- System efficiency test.
- Suitable for retinal, endovitrear, glaucoma and cataract surgery.
- Wide range of cryoprobes.



## TECHNICAL SHEET

### TECHNICAL DATA

- **GAS TYPE**  
N<sub>2</sub>O or CO<sub>2</sub> medical grade in non syphon cylinder
- **WORKING PRESSURE**  
50/60 bar (725-870psi)
- **MAX PRESSURE**  
70 bar (1015psi)
- **POWER SUPPLY**
- **INPUT VOLTAGE**  
100/120/220/240V (User Selectable)
- **FREQUENCY**  
50/60Hz
- **POWER CONSUMPTION**  
80VA

### DIMENSIONS

- **HEIGHT**  
172mm
- **WIDTH**  
393mm

- **DEPTH**  
463mm
- **WEIGHT**  
10Kg (22 lbs)
- **CLASSIFICATION AND SAFETY STANDARDS**
- **COMPLIANCE**  
EN60601-1, EN60601-1-1, EN60601-1-2, 93/42/CEE
- **CLASSIFICATION**  
Class I, Tipe BF (Applied part)
- **PROTECTION AGAINST THE INGRESS OF WATER**  
Console IPX1  
Footswitch: IPX8

### CRYO PROBES



133001 Cryo probe for cataract



133002 Cryo probe for retina, ball



133003 Cryo probe for retina, spatula



133004 Cryo probe for retina, "T" shaped



133005 Cryo probe for retina, "Bonnet"



133006 Cryo probe for glaucoma



133007 Cryo probe, endovitrear

### MISCELLANEOUS



132005 Spare cartridge for gas filter



132006 CO<sub>2</sub> cylinder adapter (Chinese standard)



132008 N<sub>2</sub>O cylinder adapter (ISO5145 standard)



132009 CO<sub>2</sub> cylinder adapter (ISO5145 standard)



132010 N<sub>2</sub>O cylinder adapter (French standard)



183001 Trolley for Cryoline with cylinder holder

■ **OPTIKON 2000 S.p.a.**

Via del Casale di Settebagni, 13  
00138 Rome  
[www.optikon.com](http://www.optikon.com)

■ **CONTACTS**

■ **CUSTOMER SERVICE**

**tel.** +39. 06.88.88.355  
**fax** +39. 06.88.88.440  
**e-mail:** [sales@optikon.com](mailto:sales@optikon.com)

■ **WORLD DEALERS**

**EUROPE**

Austria  
Belgium  
Denmark  
France  
Germany  
Greece  
Italy  
Netherlands  
Poland  
Portugal  
United Kingdom  
Czech Republic  
Romania  
Russia  
Slovakia  
Spain  
Switzerland  
Turkey  
Ukraine  
Hungary

**AFRICA**

Algeria  
Egypt  
Morocco  
South Africa  
Sudan  
Tunisia

**NORTH AMERICA**

Canada  
USA

**CENTRAL AMERICA**

Cuba  
Guatemala  
Mexico

**SOUTH AMERICA**

Argentina  
Brasil  
Chile  
Colombia  
Ecuador  
Peru

**MIDDLE EAST**

Saudi Arabia  
Bahrain  
Jordan  
Iran  
Iraq  
Lebanon  
Pakistan  
Syria

**ASIA**

Bangladesh  
China  
South korea  
Japan  
Hong kong  
India  
Indonesia  
Nepal  
Philippines  
Thailand  
Taiwan  
Vietnam



Distributed by:

TM denotes trademark of Optikon 2000 S.p.A.  
All right reserved.

