



20G Biconical Phaco tip		Sleeve
INC*	15°	
2.8	113230	113308



21G Biconical Phaco tip		Sleeve
INC*	15°	
2.3	113233	113307



Standard Phaco tip 21G**		Sleeve		
INC*	5°	15°	30°	
2.3	113235	113236	113237	113307



19G Bordin "High Followability" Phaco tip		Sleeve
INC*		
3.2	113232	113301



20G Bordin "High Followability" Phaco tip		Sleeve
INC*		
2.8	113229	113308



22G Carbon coat + diaphragm Phaco tip		Sleeve			
INC*	5°	30°	45°	Progressive	
1.0	113208	113209	113210	113211	No sleeve: bimanual Phaco



113217 Pars plana Phaco tip

* Recommended Incision Size: mm

SLEEVES



113301 Silicone sleeve for 19G U/S tips (3pcs), incision size 3.2 mm



113308 Silicone sleeve for 20G U/S tips (3pcs), incision size 2.8 mm



113307 Silicone sleeve for 21G U/S tips (3pcs), incision size 2.3 mm



113303 PTFE sleeve for 19G U/S tips, incision size 3.2 mm (3pcs)

113304 PTFE sleeve for 20G U/S incision size 2.8 mm

ACCESSORIES



116001 Test chamber (2pcs)



113401 Phaco tip wrench



113402 Phaco tip wrench



191004 Kiss (Kawano Irrigation Surge Suppressor), single-use (20 pcs)



Pulsar



OPTIKON 2000 S.p.a.

Via del Casale di Settebagni, 13
00138 Rome
www.optikon.com

CONTACTS

CUSTOMER SERVICE

tel. +39. 06.88.88.355
fax +39. 06.88.88.440
e-mail: sales@optikon.com

WORLD DEALERS

EUROPE

Austria
Belgium
Denmark
France
Germany
Greece
Italy
Netherlands
Poland
Portugal
United Kingdom
Czech Republic
Romania
Russia
Slovakia
Spain
Switzerland
Turkey
Ukraine
Hungary

AFRICA

Algeria
Egypt
Morocco
South Africa
Sudan
Tunisia

NORTH AMERICA

Canada
USA

CENTRAL AMERICA

Cuba
Guatemala
Mexico

SOUTH AMERICA

Argentina
Brasil
Chile
Colombia
Ecuador
Peru

MIDDLE EAST

Saudi Arabia
Bahrain
Jordan
Iran
Iraq
Lebanon
Pakistan
Syria

ASIA

Bangladesh
China
South Korea
Japan
Hong Kong
India
Indonesia
Nepal
Philippines
Thailand
Taiwan
Vietnam



A11 1001/01/10/English



Distributed by:

TM denotes trademark of Optikon 2000 S.p.a.
All right reserved.

Pulsar
Minimal Stress

ACCESSORIES



FORCEPS PROBES HANDPIECES TIPS SLEEVES ACCESSORIES

FORCEPS



116006 Capsulorhexis forceps

PROBES



114303 Irrigating nucleus manipulator



114304 Irrigating chopper

HANDPIECES



113101 Phaco handpiece "Slim 4"



Phaco handpiece single-use (2 crystals, single-use, sterile, 5 pcs, requires adaptor REF 116005)

RIS*	5°	30°	45°	Diaphragm tip	Progressive tip
3.2	191311	191313	191314	191315	191316

116005 Adapter for single-use handpiece



TIPS



19G Standard Phaco tip

INC*	5°	15°	30°	45°	Sleeve
3.2	113201	113202	113203	113204	113301-113303



19G Kawano double claw Phaco tip

INC*	15°	30°	Hybrid	Sleeve
3.2	113212	113213	113214	113301-113303



19G Pezzola diaphragm Phaco tip

INC*	30°	Sleeve
3.2	113218	113301-113303



20G Standard Phaco tip

INC*	5°	15°	30°	45°	Sleeve
2.8	113234	113205	113206	113207	113308-113304



20G Carbon coat + diaphragm Phaco tip

INC*	5°	15°	30°	Angled	Sleeve
1.5	113219	113220	113221	113222	No sleeve: bimanual phaco



19G Biconical Phaco tip

INC*	15°	Sleeve
3.2	113228	113301

PROBES CANNULAS HANDPIECES VITREOUS CUTTERS/MICROSCISSORS

PROBES



114301 Irrigation probe for bimanual technique



114302 Aspiration probe for bimanual technique

CANNULAS



19G I/A cannula with silicone sleeve

	INC*	
Straight	3.2	114201
Angled	3.2	114203



20G I/A cannula with silicone sleeve

	INC*	
Straight	2.8	114216
Angled	2.8	114217



21G I/A cannula with silicone sleeve

	INC*	
Straight	2.3	114218
Angled	2.3	114219



I/A cannula with metal sleeve

	INC*	
Straight	3.2	114205
Angled	3.2	114207



114211 12 o'clock I/A cannula, metal sleeve, incision size 3.0 mm



114212 90° I/A cannula, metal sleeve, incision size 3.0 mm

HANDPIECES



114101 I/A handpiece, long barrel



114102 I/A handpiece, short barrel

VITREOUS CUTTERS/MICROSCISSORS



114401 Sleeve for anterior vitrectomy (3pcs)



123002 Anterior vitrectomy probe, w/silicone sleeve

On demand single-use phaco packs for Pulsar (10pcs) including I/A tubing set for Pulsar, screen drape, Mayo drape, test chamber, sleeve (to customer choice).



* Recommended Incision Size: mm

TUBING DIATHERMY MISCELLANEA

TUBING



116002 Drainage bag (10pcs)



117001 I/A tubing set with ACS3 device for Pulsar



117002 ACS3 device for Pulsar



191701 I/A tubing set for Pulsar (single-use, sterile, 10pcs, to be used with ACS3 device REF. 117002)



191702 I/A tubing set for Pulsar with ACS3 device and drainage bag (single-use, sterile, 10pcs)

DIATHERMY



115101 Diathermy bipolar cable



115301 Diathermy forceps



115302 Diathermy pencil eraser

MISCELLANEA



112101 Dual linear footpedal




181003 I.V. pole extension including blocking screw



181004 Equipped motorized cart



191602 Screen & Mayo drape single-use (5pcs)



116003 Sterilization tray 25x15x38 cm

* Recommended Incision Size: mm