

Pulsar 2

Minimal Stress
& VitModule



- ✓ **Double pump (Venturi & Peristaltic)**
- ✓ **Minimal IOP forced irrigation**
- ✓ **Minimal Stress for M.I.C.S. & co- M.I.C.S.**
- ✓ **23G & 25G High speed vitrectomy**
- ✓ **HID Xenon light**
- ✓ **Programmable double linear footswitch**

customer support (+39)06 88.88.355
fax (+39)06 88.88.440

e-mail sales@optikon.com
www.optikon.com



OPTIKON™
OPHTHALMIC EQUIPMENT